

**Technical note T154:** Different materials step height – new method

# Error correction in different material step height measurement using Coherence Correlation Interferometry

Yang Yu, PhD Senior Applications Scientist

This Technical note shows how CCI can resolve issues with ‘different material step height measurement errors’ explained in Technical note ‘T153’.

Both thin film steps and dissimilar material steps give rise to a phase-change on reflection (PCOR). This inevitably distorts the interference series and result in the presence of a DC-shift in the measured topographical surface, leading to an error in step height measurement. CCI will resolve these issues.

## 1 PCOR induced error correction – ‘Films and Materials’ technique

### ‘Films and Materials’ technique

Through the knowledge of thin film structure and ‘Films and Materials’ technique, the DC shift (shift of the fringe envelope) due to the phase-change on reflection (PCOR) for both dissimilar material step and thin film step can be compensated so that the ‘true’ step-heights are determined.

### Error correction

A pair of HCF [1,2] functions are directly determined from the CCI interference series corresponding to the top surface and the bottom surface of the step.

A corresponding pair of synthetic functions is then optimised to match this function pair. These synthetic functions take account of the spectral n and k values of the substrate (bottom of step) and likewise of any thin film or dissimilar material present (top of step); Accurate step height is then possible through taking advantage of the topographical terms embedded within these synthetic functions.

## 2 Measurement range and accuracy

### Accuracy – ~2%

‘Films and Materials’ technique can be used for the step height measurement of dissimilar material step and thin film step with an accuracy of ~2% accuracy

### Measurement range

Dissimilar material step	0 to ~130 um (latest software)
Thin film step	~25 nm to 5 µm

The measurement can be made within the measurement range of ~25 nm to 5 µm for a step with a semi-transparent thin film coated surface, and ~ 0 to scan range limited step height value for a dissimilar material step.

### 3 How it works on CCI HD?



Figure 2(a) An example of a single layer thin film step height measurement use CCI HD

2. Top surface of the step

1. Bottom surface of the step

Region:	Name:	Bare Substrate:	Material / Film ID:	Ref:	Step Index:	Diag:	Left:	Rectangle Top:	Width:	Height:
1	Bottom	<input checked="" type="checkbox"/>	B270	<input type="checkbox"/>	1	<input checked="" type="checkbox"/>	326	472	38	40
2	Top	<input type="checkbox"/>		<input type="checkbox"/>	2	<input checked="" type="checkbox"/>	320	397	50	36



CCI Film Thickness Results


Action	Region	Thickness 1	Thickness 2	Thickness 3	Step Height	Fit	X	Y	Status	Date
Reference	ref	0.0	0.0	0.0	0.0	0.000	68.777	47.870	OK	19-Jun-2013 12:26:51
Film	Bottom	0.0	0.0	0.0	0.0	0.007	68.777	47.870	OK	19-Jun-2013 12:26:54
Film	Top	493.5	0.0	0.0	478.2	0.501	68.777	47.870	OK	19-Jun-2013 12:27:04

Film thickness – an additional benefit is the film thickness value for each layer is also obtained while calculating the step height.

Step height result

## 4 Comparison results

The 6 reference thin films were used in this study, the dissimilar material step height results were compared with other techniques include spectrophotometry and film thickness measurement



### Step height comparison results

Thin film step (reference films)

	Film 1	Film 2	Film 3	Film 4	Film 5	Film 6
Film/substrate	Ta205/ B270	Ta205/ B270	SiO <sub>2</sub> /Si	Ta205/ B270	Ta205/ B270	SiO <sub>2</sub> / Si
Film thickness (spectrophotometer) (nm)	48.8	98.5	311.1	504.6	1002.2	994.5
Mean film thickness (film thickness measurement) (nm)	52.5 ± 0.6	99.1 ± 0.3	312.8 ± 0.3	505.1 ± 0.6	1003.4 ± 1.2	996.9 ± 1.2
Mean step height (dissimilar material step measurement) (nm)	<b>56.1 ± 0.7</b>	<b>97.3 ± 0.3</b>	<b>315.9 ± 2.1</b>	<b>514.1 ± 1.1</b>	<b>994.8 ± 1.6</b>	<b>1014.6 ± 1.6</b>

The standard deviation ( $\sigma$ ) for all the analysis was calculated by means of 20 measurements


**Good agreement!**

## 4 Summary

Through the development of 'Films and Materials' technique, coherence correlation interferometry (CCI) has become the ideal method to obtain not only film thicknesses and interface information but also true step heights for dissimilar material steps and thin film steps.

'Films and Materials' technique can eliminate the step height measurement errors caused by the distortion of the fringe series due to PCOR from both thin film steps and dissimilar material steps.

## Following up technical note T153




## Some example applications



Mems devices



Solar PV



Semiconductors



Hybrid circuits and electrical contacts

## References

- 1 Mansfield, D., 'The distorted helix: thin film extraction from scanning white light Interferometry,' Proc. SPIE vol6186, 2006
- 2 Mansfield, D., US7755768, 'Apparatus for and a method of determining a characteristic of a layer or layers', 2010



**Taylor Hobson UK**  
(Global Headquarters)  
PO Box 36, 2 New Star Road  
Leicester, LE4 9JD, England  
Tel: +44 116 276 3771  
taylor-hobson.sales@ametek.com



**Taylor Hobson France**  
Tel: +33 130 68 89 30  
taylor-hobson.france@ametek.com



**Taylor Hobson Germany**  
Tel: +49 611 973040  
taylor-hobson.germany@ametek.com



**Taylor Hobson India**  
Tel: +91 80 67823200  
taylor-hobson.india@ametek.com



**Taylor Hobson Italy**  
Tel: +39 02 946 93401  
taylor-hobson.italy@ametek.com



**Taylor Hobson Japan**  
Tel: +81 36809 2406  
taylor-hobson.japan@ametek.com



**Taylor Hobson Korea**  
Tel: +82 31 888 5255  
taylor-hobson.korea@ametek.com



**Taylor Hobson China Beijing Office**  
Tel: +86 10 8526 2111  
taylor-hobson.beijing@ametek.com



**Taylor Hobson China Shanghai Office**  
Tel: +86 21 58685111-110  
taylor-hobson.shanghai@ametek.com



**Taylor Hobson Singapore**  
Tel: +65 6484 2388 Ext 120  
taylor-hobson.singapore@ametek.com



**Taylor Hobson USA**  
Tel: +1 630 621 3099  
taylor-hobson.usa@ametek.com